

## Human CD16b (FcγR3b) protein (Recombinant) (STJP000294)

STJP000294

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Human CD16b (FcγR3b)-protein was developed from hek293. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant CD16b is supplied as a lyophilized 0.2 μg/ml filtered PBS solution, pH7.4 containing 5% trehalose.
<b>Purification</b>	
<b>Dilution Range</b>	>95%, as determined by SDS-PAGE and HPLC
<b>Storage</b>	Recombinant CD16b, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
<b>Instruction</b>	Avoid repeated freeze/thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	2215
<b>Gene Symbol</b>	FCGR3B
<b>Uniprot ID</b>	FCG3B_HUMAN
<b>Immunogen</b>	
<b>Sequence</b>	