

Human alpha 1 antitrypsin protein (Recombinant) (STJP000304)

STJP000304

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Proteins |
| Short Description | Recombinant-Human alpha 1 antitrypsin-protein was developed from hek293. For use in research applications. |
| Host/Source | HEK293 |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Concentration | |
| Formulation | Recombinant Human alpha 1 antitrypsin is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose. |
| Purification | |
| Dilution Range | >95%, as determined by SDS-PAGE and HPLC |
| Storage Instruction | Recombinant Serpin A1 can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|------------|
| Gene ID | 5265 |
| Gene Symbol | SERPINA1 |
| Uniprot ID | A1AT_HUMAN |
| Immunogen Sequence | |