

Human alpha 1 antichymotrypsin protein (Recombinant) (STJP000305)

STJP000305

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human alpha 1 antichymotrypsin-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant human alpha 1 antichymotrypsin (AACT) is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant Human alpha 1 antichymotrypsin (AACT) can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C
Instruction	for twelve months. Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	12
Gene Symbol	SERPINA3
Uniprot ID	AACT_HUMAN
Immunogen	
Sequence	