

Human Cystatin A protein (Recombinant) (STJP000307)

STJP000307

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Cystatin A-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Cystatin A is lyophilized from 0.2 M filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant Human Cystatin A (CSTA) can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1475
Gene Symbol	CSTA
Uniprot ID	CYTA_HUMAN
Immunogen	
Sequence	