

Human TGF Beta-1 protein (Recombinant) (Fc-Tag) (STJP000315)

STJP000315

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human TGF Beta-1-Fc-Tag protein was developed from hek293 and has a target region of Fc-Tag. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant TGF-b1 is lyophilized from 0.2 Mu m filtered 100mM glycine PH 3.0.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC.
Storage	Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at-20°C for six months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7040
Gene Symbol	TGFB1
Uniprot ID	TGFB1_HUMAN
Immunogen	
Sequence	