

Human Interferon-beta 1a protein (Recombinant) (STJP000323)

STJP000323

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Interferon-beta 1a-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant IFN-bawas lyophilized from a.2 Mu m filtered sodium acetate solution pH.8.
Purification	
Dilution Range	>98%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least years from date of receipt at-20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at-20°C for six months, with a carrier protein without detectable loss
Instruction	

TARGET INFORMATION

Gene ID	3456
Gene Symbol	IFNB1
Uniprot ID	IFNB_HUMAN
Immunogen	
Sequence	