

Human IL2RA (CD25) protein (Recombinant) (STJP000332)

STJP000332

GENERAL INFORMATION


Product Type	Proteins
Short Description	Recombinant-Human IL2RA (CD25)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant CD25 is supplied as a 0.2 μ m filtered PBS solution, pH7.4.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	Recombinant CD25, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
Instruction	Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	3559
Gene Symbol	IL2RA
Uniprot ID	IL2RA_HUMAN
Immunogen Sequence	



sIL-2R Alpha treated mice displayed increased novelty-induced stereotypic behaviors. In particular, sIL-2R Alpha treated mice showed dose-dependent increased head-bobbing