

Human Serum Albumin protein (Recombinant) (STJP000336)

STJP000336

GENERAL INFORMATION

Product Type Proteins
Short Description Recombinant-Human Serum Albumin-protein was developed from p.pastoris. For use in research applications.
Host/Source P.pastoris

PRODUCT PROPERTIES

Concentration
Formulation Recombinant HSA was lyophilized from mg/ml. 2 µm filtered solution in PBS, pH.
Purification
Dilution Range >99%, as determined by SDS-PAGE and HPLC
Storage The lyophilized protein is stable for at least years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss
Instruction

TARGET INFORMATION

Gene ID 213
Gene Symbol ALB
Uniprot ID ALBU_HUMAN
Immunogen
Sequence