

HIV-1 gp41 protein (Recombinant) (Full Length) (STJP000345)

STJP000345

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-HIV-1 gp41-Full Length protein was developed from e.coli and has a target region of Full Length. For use in research applications.
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	
Formulation	20mM PBS pH-7.8, 20mM NaCl 0.5M, 1mM DTT & 8M urea.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE.
Storage	This protein can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	