

Human GM-CSF protein (Recombinant) (STJP000349)

STJP000349

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human GM-CSF-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant CSF2 was lyophilized from 0.2 µm filtered PBS solution pH7.4
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage Instruction	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C.

TARGET INFORMATION

Gene ID	1437
Gene Symbol	CSF2
Uniprot ID	CSF2_HUMAN
Immunogen Sequence	