

## Mouse TGFB1 protein (Recombinant) (STJP000350)

STJP000350

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Mouse TGFB1-protein was developed from cho cells. For use in research applications.
<b>Host/Source</b>	CHO cells

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant TGFB1 was lyophilized from 0.2 µm filtered PBS solution pH7.4
<b>Purification</b>	
<b>Dilution Range</b>	>95%, as determined by SDS-PAGE and HPLC
<b>Storage Instruction</b>	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C.

### TARGET INFORMATION

<b>Gene ID</b>	21803
<b>Gene Symbol</b>	Tgfb1
<b>Uniprot ID</b>	TGFB1_MOUSE
<b>Immunogen Sequence</b>	