

Mouse Interleukin-5 protein (Recombinant) (STJP000354)

STJP000354

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse Interleukin-5-protein was developed from insect cells. For use in research applications.
Host/Source	Insect cells

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Interleukin-5 was lyophilized from a.2 Mu m filtered PBS solution.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss
Instruction	

TARGET INFORMATION

Gene ID	16191
Gene Symbol	IL5
Uniprot ID	IL5_MOUSE
Immunogen	
Sequence	