

Human HIV-1 gp120 protein (Recombinant) (STJP000355)

STJP000355

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human HIV-1 gp120-protein is for use in research applications.
Host/Source	

PRODUCT PROPERTIES

Concentration	
Formulation	0.2 µm filtered 50mM Tris pH 8.0 NaCl 60mM, 10mM Glutathione, 0.25% Sarkosyl, 50% Glycerol
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	This product can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	