

Mouse Soluble CD4 protein (Recombinant) (STJP000364)

STJP000364

GENERAL INFORMATION

Product Type Proteins
Short Description Recombinant-Mouse Soluble CD4-protein was developed from hek293. For use in research applications.
Host/Source HEK293

PRODUCT PROPERTIES

Concentration
Formulation Recombinant CD4 was lyophilized from a.2 Mu m filtered solution inPBS, pH.4.
Purification
Dilution Range >90%, as determined by SDS-PAGE and HPLC
Storage The lyophilized protein is stable for at least years from date of receipt at-20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at-20°C for six months, with a carrier protein without detectable loss
Instruction

TARGET INFORMATION

Gene ID [12504](#)
Gene Symbol [Cd4](#)
Uniprot ID [CD4_MOUSE](#)
Immunogen
Sequence