

Human CD66c (CEACAM6) protein (Recombinant) (STJP000369)

STJP000369

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CD66c (CEACAM6)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant human CD66cis supplied as a 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant CD66c, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
Instruction	Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	1048
Gene Symbol	CEACAM5
Uniprot ID	CEAM5_HUMAN
Immunogen Sequence	