

Human CSF-1 receptor protein (Recombinant) (STJP000375)

STJP000375

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CSF-1 receptor-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant CSF-1R is supplied as lyophilized 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>98%, as determined by SDS-PAGE and HPLC
Storage	Recombinant CSF-1R, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1436
Gene Symbol	CSF1R
Uniprot ID	CSF1R_HUMAN
Immunogen Sequence	