

Human Cathepsin B protein (Recombinant) (STJP000382)

STJP000382

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Cathepsin B-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage Instruction	

TARGET INFORMATION

Gene ID	1508
Gene Symbol	CTSB
Uniprot ID	CATB_HUMAN
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081