

Human LYN protein (Recombinant) (STJP000383)

STJP000383

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human LYN-protein was developed from insect cells. For use in research applications.
Host/Source	Insect cells

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human LYN was lyophilized from a 0.2 µg/ml filtered 20mM Tris, 50mM NaCl, pH 8.0, 0.5mM PMSF.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE
Storage	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss.
Instruction	

TARGET INFORMATION

Gene ID	4067
Gene Symbol	LYN
Uniprot ID	LYN_HUMAN
Immunogen	
Sequence	