

Human Activin protein (Recombinant) (STJP000387)

STJP000387

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Activin-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	
Formulation	Human Activin A was lyophilized from a 0.2 Mu m filtered solution in 0.2M AcOH, pH 4.0.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable los
Instruction	

TARGET INFORMATION

Gene ID	3624
Gene Symbol	INHBA
Uniprot ID	INHBA_HUMAN
Immunogen Sequence	