

Rat CD2 (LFA-2) protein (Recombinant) (STJP000394)

STJP000394

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Rat CD2 (LFA-2)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Rat CD2 is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	Recombinant Rat CD2 can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	497761
Gene Symbol	Cd2
Uniprot ID	CD2_RAT
Immunogen	
Sequence	