

Mouse CSF-1 receptor (CSF-1R) protein (Recombinant) (STJP000402)

STJP000402

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse CSF-1 receptor (CSF-1R)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant mouse CSF-1R is supplied as lyophilized 0.2 µg/ml filtered PBS solution, pH7.2.
Purification	
Dilution Range	>98%, as determined by SDS-PAGE and HPLC
Storage	Recombinant CSF-1R, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	12978
Gene Symbol	Csf1r
Uniprot ID	CSF1R_MOUSE
Immunogen	
Sequence	