

Human PDGFR beta protein (Recombinant) (STJP000403)

STJP000403

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human PDGFR beta-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant PDGFR beta is supplied as a 0.2 μ g/ml filtered PBS solution, pH7.2.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant PDGFRb, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5159
Gene Symbol	PDGFRB
Uniprot ID	PGFRB_HUMAN
Immunogen	
Sequence	