

Human Placental lactogen protein (Recombinant) (STJP000407)

STJP000407

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Placental lactogen-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant human Placental lactogen (CSH2) is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant Human Placental lactogen (CSH1) can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	51060
Gene Symbol	TXNDC12
Uniprot ID	TXD12_HUMAN
Immunogen	
Sequence	