

Mouse Lipocalin-2 protein (Recombinant) (STJP000411)

STJP000411

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse Lipocalin-2-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Mouse Lipocalin-2 is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant NGAL can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	16819
Gene Symbol	Lcn2
Uniprot ID	NGAL_MOUSE
Immunogen	
Sequence	