

## Human CD32a (FcγRIIa) protein (Recombinant) (STJP000413)

STJP000413

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Human CD32a (FcγRIIa)-protein was developed from hek293. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant CD32a is supplied as a lyophilized 0.2 μg/ml filtered PBS solution, pH7.2.
<b>Purification</b>	
<b>Dilution Range</b>	>95%, as determined by SDS-PAGE and HPLC
<b>Storage</b>	Lyophilized CD32a, as supplied, should be stored at -20°C to -70°C. Upon reconstitution, working aliquots may be stored at 2-8°C for one month, or six months at -20°C without detectable loss of activity. It is recommended to add 0.1% BSA for long term storage.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	2212
<b>Gene Symbol</b>	FCGR2A
<b>Uniprot ID</b>	FCG2A_HUMAN
<b>Immunogen</b>	
<b>Sequence</b>	