

## Human CD66a (CEACAM1) protein (Recombinant) (STJP000421)

STJP000421

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Human CD66a (CEACAM1)-protein was developed from hek293. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant CD66ais supplied as a 0.2 Mu m filtered PBS solution, pH7.2.
<b>Purification</b>	
<b>Dilution Range</b>	>95%, as determined by SDS-PAGE and HPLC
<b>Storage</b>	Recombinant CD66a, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
<b>Instruction</b>	Avoid repeated freeze/thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">634</a>
<b>Gene Symbol</b>	<a href="#">CEACAM1</a>
<b>Uniprot ID</b>	<a href="#">CEAM1_HUMAN</a>
<b>Immunogen Sequence</b>	