

mouse EPO receptor (EPOR) protein (Recombinant) (STJP000429)

STJP000429

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-mouse EPO receptor (EPOR)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant mouse EPO receptor is supplied as a 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	Recombinant mouse EPOR, as supplied, can be stored in working aliquots at 2-8°C for one month, or at-20°C to-70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	13857
Gene Symbol	Epor
Uniprot ID	EPOR_MOUSE
Immunogen	
Sequence	