

Human IL-1 receptor 1 (CD121a) protein (Recombinant) (STJP000430)

STJP000430

GENERAL INFORMATION


Product Type Proteins
Short Description Recombinant-Human IL-1 receptor 1 (CD121a)-protein was developed from hek293. For use in research applications.
Host/Source HEK293

PRODUCT PROPERTIES

Concentration
Formulation Recombinant Human IL-1R (CD121a) is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification
Dilution Range >97%, as determined by SDS-PAGE and HPLC
Storage Instruction Recombinant Human IL1R1 can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID [3554](#)
Gene Symbol [IL1R1](#)
Uniprot ID [IL1R1_HUMAN](#)
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081