

Human Carboxypeptidase A1 protein (Recombinant) (STJP000433)

STJP000433

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Carboxypeptidase A1-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Carboxypeptidase A1 is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant Carboxypeptidase A1 (CPA1) can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1357
Gene Symbol	CPA1
Uniprot ID	CBPA1_HUMAN
Immunogen	
Sequence	