

Human CD155 protein (Recombinant) (Fc-Tag) (STJP000434)

STJP000434

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CD155-Fc-Tag protein was developed from hek293 and has a target region of Fc-Tag. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant CD155is supplied as a 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant CD155, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
Instruction	Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	5817
Gene Symbol	PVR
Uniprot ID	PVR_HUMAN
Immunogen Sequence	