

Human CD119 (IFNGR1) protein (Recombinant) (STJP000438)

STJP000438

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CD119 (IFNGR1)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human IFNGR1 (CD119) is supplied as a 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage Instruction	Recombinant Human IFN gamma receptor (CD119) , as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	3459
Gene Symbol	IFNGR1
Uniprot ID	INGR1_HUMAN
Immunogen Sequence	