

Human TIMP-2 protein (Recombinant) (STJP000445)

STJP000445

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human TIMP-2-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human TIMP-2 is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant TIMP2 can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7077
Gene Symbol	TIMP2
Uniprot ID	TIMP2_HUMAN
Immunogen	
Sequence	