

Human CD31 (PECAM1) protein (Recombinant) (Fc-Tag) (STJP000449)

STJP000449

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CD31 (PECAM1)-Fc-Tag protein was developed from hek293 and has a target region of Fc-Tag. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant human CD31 is supplied as a 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant CD31, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
Instruction	Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	5175
Gene Symbol	PECAM1
Uniprot ID	PECA1_HUMAN
Immunogen Sequence	