

Human CD36 (FAT) protein (Recombinant) (STJP000453)

STJP000453

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CD36 (FAT)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human CD36 is lyophilized from 0.2 M filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>90%, as determined by SDS-PAGE and HPLC
Storage	Recombinant CD36 can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	948
Gene Symbol	CD36
Uniprot ID	CD36_HUMAN
Immunogen	
Sequence	