

Mouse Cathepsin D protein (Recombinant) (STJP000461)

STJP000461

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse Cathepsin D-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant mouse Cathepsin D is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage Instruction	Recombinant Cathepsin D can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	13033
Gene Symbol	Ctsd
Uniprot ID	CATD_MOUSE
Immunogen Sequence	