

## Mouse L-Selectin (CD62L) protein (Recombinant) (STJP000462)

STJP000462

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Mouse L-Selectin (CD62L)-protein was developed from hek293. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant Mouse L-Selectin (CD62L) is lyophilized from 0.2 M filtered PBS solution, pH7.2, 5% Trehalose.
<b>Purification</b>	
<b>Dilution Range</b>	>97%, as determined by SDS-PAGE and HPLC
<b>Storage Instruction</b>	Recombinant L-Selectin can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	20343
<b>Gene Symbol</b>	Sell
<b>Uniprot ID</b>	LYAM1_MOUSE
<b>Immunogen</b>	
<b>Sequence</b>	