

Human CD106 (VCAM1) protein (Recombinant) (STJP000467)

STJP000467

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CD106 (VCAM1)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant CD106is supplied as a 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant CD106, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
Instruction	Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	7412
Gene Symbol	VCAM1
Uniprot ID	VCAM1_HUMAN
Immunogen Sequence	