

Human Follistatin protein (Recombinant) (STJP000469)

STJP000469

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Follistatin-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human Follistatin was lyophilized at 1mg/ml containing no additives.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss.
Instruction	

TARGET INFORMATION

Gene ID	10468
Gene Symbol	FST
Uniprot ID	FST_HUMAN
Immunogen	
Sequence	