

## Human Transcobalamin II protein (Recombinant) (STJP000470)

STJP000470

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Human Transcobalamin II-protein was developed from hek293. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant human Transcobalamin II is lyophilized from 0.2 M unfiltered PBS solution, pH7.2, 5% Trehalose.
<b>Purification</b>	
<b>Dilution Range</b>	>95%, as determined by SDS-PAGE and HPLC
<b>Storage</b>	Recombinant Transcobalamin II can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
<b>Instruction</b>	Avoid repeated freeze/thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">6948</a>
<b>Gene Symbol</b>	<a href="#">TCN2</a>
<b>Uniprot ID</b>	<a href="#">TCO2_HUMAN</a>
<b>Immunogen Sequence</b>	