

Mouse TNFR1 (CD120a) protein (Recombinant) (STJP000489)

STJP000489

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse TNFR1 (CD120a)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Mouse TNFR1 (TNFRSF1A) is supplied as a 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant Mousen TNF receptor 1 (CD120a) , as supplied, can be stored in working aliquots at 2-8°C for one month, or at-20°C
Instruction	to-70°C for twelve months. Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	21937
Gene Symbol	Tnfrsf1a
Uniprot ID	TNR1A_MOUSE
Immunogen Sequence	