

## Mouse TNFR2 (CD120b) protein (Recombinant) (STJP000490)

STJP000490

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Mouse TNFR2 (CD120b)-protein was developed from hek293. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant Mouse TNFR2 (TNFRSF1B) is supplied as a 0.2 µg/ml filtered PBS solution, pH7.2.
<b>Purification</b>	
<b>Dilution Range</b>	>95%, as determined by SDS-PAGE and HPLC
<b>Storage</b>	Recombinant Mouse TNF receptor 2 (CD120b) , as supplied, can be stored in working aliquots at 2-8°C for one month, or at-20°C to-70°C for twelve months. Avoid repeated freeze/thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">21938</a>
<b>Gene Symbol</b>	<a href="#">Tnfrsf1b</a>
<b>Uniprot ID</b>	<a href="#">TNR1B_MOUSE</a>
<b>Immunogen Sequence</b>	