

Mouse TGFB2 protein (Recombinant) (STJP000495)

STJP000495

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse TGFB2-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant TGFB2 was lyophilized from 0.2 µm filtered PBS solution pH7.4
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage Instruction	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C.

TARGET INFORMATION

Gene ID	21808
Gene Symbol	Tgfb2
Uniprot ID	TGFB2_MOUSE
Immunogen Sequence	