

Mouse CD40 (TNFRSF5) protein (Recombinant) (STJP000497)

STJP000497

GENERAL INFORMATION


Product Type	Proteins
Short Description	Recombinant-Mouse CD40 (TNFRSF5)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Mouse CD40 is lyophilized from 0.2 M unfiltered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	Recombinant CD40 can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	21939
Gene Symbol	Cd40
Uniprot ID	TNRS_MOUSE
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081