

Mouse IL-1 receptor type II (CD121b) protein (Recombinant) (STJP000500)

STJP000500

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse IL-1 receptor type II (CD121b)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Mouse IL-1 receptor II (CD121b) is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose, 0.01% tween80
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	Recombinant Mouse IL1R2 can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	16178
Gene Symbol	Il1r2
Uniprot ID	IL1R2_MOUSE
Immunogen	
Sequence	