

Human CD30 (TNFRSF8) protein (Recombinant) (STJP000503)

STJP000503

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CD30 (TNFRSF8)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human CD30 is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	Recombinant CD30 can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	943
Gene Symbol	TNFRSF8
Uniprot ID	TNR8_HUMAN
Immunogen	
Sequence	