

## Rat Lipocalin-2 protein (Recombinant) (STJP000508)

STJP000508

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Rat Lipocalin-2-protein was developed from hek293. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant Rat Lipocalin-2 is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
<b>Purification</b>	
<b>Dilution Range</b>	>95%, as determined by SDS-PAGE and HPLC
<b>Storage</b>	Recombinant NGAL can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">170496</a>
<b>Gene Symbol</b>	<a href="#">Lcn2</a>
<b>Uniprot ID</b>	<a href="#">NGAL_RAT</a>
<b>Immunogen</b>	
<b>Sequence</b>	