

Human CD80 (B7-1) protein (Recombinant) (STJP000516)

STJP000516

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CD80 (B7-1)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant CD80 was lyophilized from 0.2µm filtered 50mM Tris, 100mM NaCl, pH 8.0.
Purification	
Dilution Range	>90%, as determined by SDS-PAGE and HPLC
Storage	Recombinant CD80, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
Instruction	Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	941
Gene Symbol	CD80
Uniprot ID	CD80_HUMAN
Immunogen Sequence	