

Human CD68 protein (Recombinant) (STJP000517)

STJP000517

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CD68-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant human CD68 (macrosialin) is supplied as lyophilized 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>98%, as determined by SDS-PAGE and HPLC
Storage	Recombinant human CD68 (macrosialin) , as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C
Instruction	for twelve months. Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	968
Gene Symbol	CD68
Uniprot ID	CD68_HUMAN
Immunogen	
Sequence	