

Mouse IL-2RG (CD132) protein (Recombinant) (STJP000519)

STJP000519

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse IL-2RG (CD132)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human IL2RG was lyophilized from a 0.2 Mu m filtered PBS solution pH 7.4.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE
Storage	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable los
Instruction	

TARGET INFORMATION

Gene ID	16186
Gene Symbol	Il2rg
Uniprot ID	IL2RG_MOUSE
Immunogen	
Sequence	