

Mouse VEGFR2 protein (Recombinant) (Fc-Tag) (STJP000524)

STJP000524

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse VEGFR2-Fc-Tag protein was developed from hek293 and has a target region of Fc-Tag. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant KDR was lyophilized from 0.2 µm filtered PBS solution pH7.4
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage Instruction	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C.

TARGET INFORMATION

Gene ID	16542
Gene Symbol	Kdr
Uniprot ID	VGFR2_MOUSE
Immunogen	
Sequence	